

## CONFIGURACIÓN ÓPTICA CANTO II NEW

Láser	DETECTOR	FILTRO	fluorocromo clásico
blue 488nm	FL1-A	530/30	FITC, EGFP, EYFP, Alexa 488
	FL2-A	582/15	IP, PE
	FL3-A	670LP	PerCP, PerCP-Cy5.5, PE-Cy5, PE-Cy7, 7AAD
	FL4-A	780/60	PE-Cy7
red 640	FL5-A	670/30	APC, Alexa 647
	FL6-A	780/60	APC-Cy7, APC eFluor 780, APC H7, Live/Dead fixable NIR
violet 450	FL7-A	450/40	DAPI/ , Brilliant Violet 421, horizon V450
	FL8-A	605/10	Brilliant Violet 605, Ghost Dye violet 540